ACTIVE
LUXE

S · K · I · N

LED Light Therapy Face Mask

User Manual
Uporabniški priročnik
Korisnički priručnik

ACTIVE
LUXE
S · K · I · N

*LED
Light Therapy
Face Mask*

Product Description

The light energy rejuvenation beauty mask uses bioactive cold light, representing the latest in beauty technology. Operating at a very low temperature, it utilizes energy-efficient pure light and near-infrared light technology.


- Red light has a wavelength of 620-625nm, penetrates deep into the dermis, stimulating collagen and elastin production.
- Green light has a wavelength of 520-525nm, neutralizes the balance of water and oil ratio, soothing the skin.
- Blue light has a wavelength of 455-460nm, refreshes, activates, and tightens the skin.
- Yellow light is a dual frequency of red and green light; it improves oxygen exchange and brightens the skin tone.
- Purple light is a dual frequency of red and blue light, delivering the benefits of both.
- Cyan light gradually enhances cell energy and accelerates metabolism.
- Orange light helps to improve skin repair and pigmentation issues.
- Infrared red light has a wavelength of 850nm. Invisible to the eye, its wavelength is longer than red light and it can reach the deep skin.

Product Parameters

Product name	LED Light Therapy Face Mask
Product model	T70
Product size	About 321 × 215 × 5 mm
Charging specification	DC 5V / 1.0A
Power	About 3.7 W
Rated voltage	3.1V
Battery capacity	2200 mAh

Charging time	About 3 hours
Net weight	355 g

Product Structure


Key Functions

- To turn on the device, long press power key. The red light will be activated by default.
- Short press the power key to switch the nursing time; 5 minutes, 10 minutes, and 15 minutes; select as needed.
- When the device is on, short press the LED key to switch the light colors to red, green, blue, yellow, purple, cyan, and orange; select as needed.
- Short press the NIR key to turn NIR mode on and off.


Usage

1. Cleanse and dry your face.


2. Remove the light energy rejuvenation beauty mask and eye mask from the packaging. Ensure the light-emitting side of the mask faces downward, with the flat surface of the eye mask facing up and the protruding part facing down. Install it following the direction of the arrow.

3. Attach the magic tape on both sides of the mask.

4. Align the eye mask of the device with the eyes, adjust the tightness, paste and fix it with magic tape, focusing on comfort.

5. Connect the power cord of the device to the USB connecting port, long press the power key to turn it on, and select the nursing timing and light color as needed.

6. It can be used with skin care products or facial masks. Using light for irradiation, use the device for 15-30 minutes for optimal results.

Note

- It is recommended to enjoy the light rejuvenation beauty with your eyes closed during your treatment. Avoid reading, watching TV, using electronic devices, or looking directly at the light source.
- Clean the mask after use. Do not rinse it directly with water; wipe it with a damp cloth. The elastic band can be washed with water.
- Use different photon irradiation according to individual circumstances. Do not continuously switch between different light sources for irradiation at the same time.
- As the material of the mask is silicone, when opening the bag, it may have a smell. This is a normal phenomenon, and it will dissipate quickly.
- The mask has a circuit board and LED lights inside. Do not squeeze, hit, or press the mask heavily.
- After use, when on the bed before going to sleep, please put the mask back into the box. Do not leave it on the bed.
- Do not use corrosive liquids to wipe the device, such as alcohol, gasoline, nail remover, etc., as these may cause product failure, cracking, or discoloration of accessories.
- People with facial plastic surgery or subcutaneous metal implantation should not use the device. Pregnant women should use the device under the guidance of a doctor.
- Do not use this device if you are using medical electronic instruments such as a pacemaker, artificial heart-lung machine, electrocardiograph, infusion set, etc.
- Pregnant women or those who are in a physiological or breastfeeding period should not use it.
- People with heart disease or abnormal blood pressure should not use it.
- People with allergies, atopic dermatitis, skin diseases, sensitive muscles, or delicate skin, should not use it.
- If you have been under hospital treatment recently or have a past medical history, please consult your doctor before use.
- This product is intended for facial and neck care only. Do not use it for other body parts or purposes.

Common Problems and Solutions

This chapter provides an overview of the problems that may be encountered when using this product and proposes solutions. If you are unable to solve any problems with the following information, please contact the customer service center in your country or region.

Questions	Possible cause	Solutions
Unable to power on	Low battery.	Try again after fully charging.
	USB power supply is not connected properly.	Check if the power supply is plugged in properly.
LED lights cannot be adjusted	USB interface is not fully connected.	Reconnect the USB interface.
	Product damage.	Please contact after-sales service.
Skin feels dry after use	Light energy might reduce skin moisture. This is normal.	Use with a facial mask or skincare products.
Unisatisfactory results	Insufficient sun protection before and after use.	If sun protection is insufficient, ultraviolet radiation can cause an anti-blackening phenomenon. Use adequate sun protection to avoid this.
For abnormal situations or unresolved issues, contact after-sales service in your language promptly. Keep the invoice and warranty card for future product maintenance.		

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


Product Presentation

Facial style 1


Face: 60 sets of luminous matrix x 240 beam light source (4 independent luminous chips built in each matrix chip)

Facial style 2


Face: 60 sets of luminous matrix x 240 beam light source (4 independent luminous chips built in each matrix chip)

Neck and Chest

Neck and Chest: 60 sets of luminous matrix x 240 beam light source (4 independent luminous chips built in each matrix chip)

Controller


Controller: 2 USB connections, Can handle both face and neck.

Operation explanation:

- Power-on button: Hold the button for 1-2 seconds to power on/off.
- Infrared control key: Turns the infrared light on/off. Can be used with any light mode.
- E Key (energy adjustment): Short press to switch between 3 energy output levels.
- M Key (light source selection): Short press to cycle through 4 different light modes.
- T + Key (time increase): Press to extend the lighting duration.
- T - Key (time decrease): Press to shorten the lighting duration.
- Available timer settings: 10 / 15 / 20 / 25 / 30 minutes (5 levels in total).

Photon Interpretation

Red light (630 nm): Anti-aging effect; promotes collagen production, reduces skin laxity, and improves elasticity.

Blue Light (460 nm): Antibacterial action against acne-causing bacteria; reduces inflammation and scarring, increases skin oxygenation, and enhances skin brightness.

Yellow Light (590 nm): Balances skin tone and texture; stimulates red blood cell activity to reduce redness and fine lines, helps fade acne marks, and supports skin stability.

Infrared Light (850 nm): Aids skin repair; accelerates metabolism, improves blood circulation, and increases skin hydration.

Mode are: Red light, 1,2,3

Blue light, 1,2,3 File


Yellow light, 1,2,3 File

Infrared light is independently controlled and can also be used with any kind of light.

Infrared Light: In each of the three modes (red, blue, or yellow light), the infrared light can be switched on or off using the NIR button. When the controller screen displays the NIR indicator, the infrared light is active.

Installation

1. Remove Dustproof film
2. Install Eyes protector
3. Install self-buckle strap
4. Connect power and controller


Wearing Diagram

After installing the photon device, adjust the strap buckle for comfort. Then connect it to the controller and switch on the photon light.

Before use, expose your face to the light for 1 minute to test for any discomfort. Begin with gear levels 1-2 to build skin tolerance, and select the energy level according to your skin's condition. Do not increase the time or frequency blindly; adjust the duration of light usage based on the actual condition of your skin. It is recommended to use this product 3-4 times per week, with each session lasting 15-25 minutes.

Charging Method

The screen indicator flashes while charging is in progress. The indicator turns off when charging is complete.

Charge the controller when it is not connected to the Photon device.

Matters Need Attention:


1. Use the product strictly according to the instructions. Do not operate it outside of the manual, as this may cause malfunction or accidents.
2. Do not use this product for purposes other than those indicated.
3. Pregnant or breastfeeding women, as well as children, must not use this product.
4. Users with physical, sensory, or mental disorders, or those undergoing treatment for dermatological conditions (such as lupus erythematosus, photosensitive eczema, or albinism), should not use this product unless under the supervision of a responsible specialist (e.g. a dermatologist).
5. Do not use this product if you have any light-sensitive disorders. Photophobia or photosensitivity may cause severe skin reactions.
6. Do not use on the same day as photosensitive skincare ingredients.
7. Protect your eyes by wearing the eye mask correctly. It is recommended to keep eyes closed during use. If you have a hereditary eye disease, do not use this product.
8. If you feel any discomfort while using the product, stop use immediately.
9. This product is not waterproof. Do not immerse it in water or use it if it has been exposed to water.
10. Keep out of reach of infants and children. Do not allow them to use this product.
11. Do not use the product if any parts are damaged.
12. This device is for personal daily use only. Do not disassemble, modify, or repair it without authorization.
13. Recommended usage: 3-4 times per week, 15-25 minutes per session.

Product Parameter

LED lamp beads


Product name	LED Facial Mask Instrument		
LED lamp bead	240	External box size	34x23x4.5 cm
Time timing	10/15/20/25/30 min	Net weight of products	0.32 kg
Service temperature	0-40 °C	Exchange	0.3A 50/60 Hz
Adapter input	100 -240 V	Adapter output	DC 5V 3000 mA
Four wavelenghts of light	Red light was 630nm, blue light 460nm, yellow light 590nm, and near-infrared light 850nm		
Individual product list	Mask instrument, controller, charging cable, instructions, fixed strap		

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