
ACTIVE
LUXE
S · K · I · N


LED
Light Therapy
Sleeping Bag

User Manual
Uporabniški priročnik
Korisnički priručnik

Product Parameters


Model ID	T2000
Product size	180 × 90 × 30 cm
LED Source	5050SMD Single Chip
LED quantity	3600 pcs
Input Voltage	AC 85-265 V
Power consumption	280 W
Wavelength	660 nm : 850 nm = 1:1
Irradiance	0 inch > 100 mW/cm ²
Beam Angle	120°
EMF	0 μT
Single Package Size	101 × 48 × 17 cm
G.W	10.5 kg
Warranty	One Year

Controller Instructions

This LED Light Therapy Sleeping Bag has a total of 3600 pcs 5050 SMD LEDs.

Even though the 850nm LEDs look dim, they are working. You can use your phone camera to see the 850nm Near-Infrared light, or you can see it in a dark environment.

The LEDs distribution is: one row 660nm Red, one row 850nm Near-Infrared — in turn. The two Red Light Therapy Blankets have individual controllers, so you can control each blanket separately.

Product Details

Soft and high-transmittance transparent PVC plastic cover to protect LEDs and allow for easy cleaning.

Durable power cord outlet.

Long-lasting zipper.


The head side.

The foot side.

Connection Instructions

Package Contents

Red Light Therapy Sauna Bag × 1

Oxford Storage Bag × 1

Plug × 1


Controller × 2


Adapter × 1

How Does Red Light Therapy Work?


630nm Wavelength Benefits

The 630nm wavelength is ideal for skin concerns, which penetrate into the skin cells and sebaceous glands to rejuvenate the skin tone and texture by smoothing and evening out uneven pigmentation. Benefits: reducing fine lines and wrinkles; improvement in skin disorders such as psoriasis; hair regrowth; healing of superficial skin cancers; treatment of acne.


660nm Wavelength Benefits

The 660nm wavelength penetrates deeper than the 630nm. The 660nm reaches into the entire range of skin issues to promote healing and regeneration. Benefits: reduced training fatigue; inflammation and swelling following injury; neuropathic pain; improved bone healing; accelerated wound healing.


810nm (Near Infrared) Wavelength Benefits

The 810nm has the ability to extend through the skull into the brain, with benefits of improving healing and recovery, recovery from stroke in certain patients, recovery from traumatic brain injury, accelerated wound healing, and hair growth.


830nm (Near Infrared) Wavelength Benefits

This wavelength is not absorbed by the body; it is able to penetrate deeper through skin and tissue. Benefits of accelerated healing and reduced infection, improved bone recovery and growth.


850nm (Near Infrared) Wavelength Benefits

The 850nm wavelength penetrates deeper than the 830nm. The 850nm reaches into the entire range of skin issues to promote healing and regeneration. Benefits: reduced training fatigue; inflammation and swelling following injury; neuropathic pain; improved bone healing; accelerated wound healing.

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